



F E N D O

MedizinTechnik



Einweg Injektionsnadel „TipSafe“ mit Entsorgungssicherung / Disposable Injection Needle „TipSafe“ with protection for disposal

- Sehr scharfe Nadelspitze, speziell für vernarbte Mukosa geeignet
- Handgriff nach Anwendung aknickbar, sichert das Verbleiben der Nadel im Tubus
- Handgriff 2-stufig arretierbar, für präzises und leichtes Unterspritzen
- Die atraumatische Metallspitze sichert die Nadel gegen seitliches Durchstechen bei starker Endoskopabwinklung
- very sharp needle tip, special for scarred Mukosa
- handle after operation kinkable, this saves the remaining of the needle in the tubus
- handle with double lock function, for exact and easy inject
- the atraumatic metal tip protects the needle against lateral piercing during a strong endoscope bending

Artikelnummer Item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Nadel Ø needle	Nadellänge needle length	steril sterile	1 VE pack. unit
FE-D2323PN2305	2,3 x 2300 mm	2,8 mm	0,7 mm	5 mm	Ⓢ	10 Stk./pcs.