


Einweg Polypektomieschlinge, asymmetrisch, multifiler Draht, Olympus HF Anschluss

Disposable polypectomy snare, asymmetrical, multifilament wire, Olympus HF connection

| Artikelnummer item no. | Ø und Länge Ø and length | Arbeitskanal working channel | Ø Schlinge Ø snare | steril sterile | 1 VE pack. unit |
|---------------------------|-----------------------------|---------------------------------|-----------------------|-------------------|--------------------|
| FE-D2323DB25 | 2,3 x 2300 mm | 2,8 mm | 25 mm | Ⓢ | 10 Stk./pcs. |



Einweg Polypektomieschlinge, oval, multifiler Draht, Olympus HF Anschluss, stufenlos verstellbar

Disposable polypectomy snare, oval, multifilament wire, Olympus HF connection, infinitely variable

| Artikelnummer item no. | Ø und Länge Ø and length | Arbeitskanal working channel | Ø Schlinge Ø snare | steril sterile | 1 VE pack. unit |
|---------------------------|-----------------------------|---------------------------------|-----------------------|-------------------|--------------------|
| FE-D2423DT10 | 2,4 x 2300 mm | 2,8 mm | 10 mm | Ⓢ | 10 Stk./pcs. |
| FE-D2423DT20 | 2,4 x 2300 mm | 2,8 mm | 20 mm | Ⓢ | 10 Stk./pcs. |
| FE-D2423DT30 | 2,4 x 2300 mm | 2,8 mm | 30 mm | Ⓢ | 10 Stk./pcs. |