



Einweg Gallengangs-Prothese (Stent), gerade, mit Goldring und Einführhilfe

Disposable biliary duct stent, straight, with gold marker and insertion aid

Artikelnummer item no.	Ø und Länge Ø and length	Flapabstand length between flaps	steril sterile	1 VE pack. unit
FE-D10F74PRO	10 FR x 70 mm	40 mm	Ⓢ	2 Stk./pcs.
FE-D10F96PRO	10 FR x 90 mm	60 mm	Ⓢ	2 Stk./pcs.
FE-D10F118PRO	10 FR x 110 mm	80 mm	Ⓢ	2 Stk./pcs.
FE-D10F1310PRO	10 FR x 130 mm	100 mm	Ⓢ	2 Stk./pcs.
FE-D10F1512PRO	10 FR x 150 mm	120 mm	Ⓢ	2 Stk./pcs.



Einweg Führungskatheter mit Goldring und Spülansatz

Disposable guiding catheter with gold marker and rinsing unit

Artikelnummer item no.	Ø und Länge Ø and length	Ø für Prothesen Ø for prostheses	Farbe color	steril sterile	1 VE pack. unit
FE-FK6F220	6 FR x 2200 mm	10 FR/11,5 FR	weiß/white	Ⓢ	2 Stk./pcs.