


**Polypektomieschlinge, oval, multifiler Draht, Olympus HF Anschluss, Handgriff demontierbar, zusätzlicher Spülanschluss, autoklavierbar 134° C**

Polypectomy snare, oval, multifilament wire, Olympus HF connection, removable handle, additional rinse connection, autoclavable 134° C

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Größe Ø snare	unsteril non-sterile	1 VE pack. unit
FE-2423DT20	2,4 x 2300 mm	2,8 mm	20 mm	o	2 Stk. + Griff / 2 pcs. + handle
FE-2423DT30	2,4 x 2300 mm	2,8 mm	30 mm	o	2 Stk. + Griff / 2 pcs. + handle



**Polypektomieschlinge, halbmond, multifiler Draht, Olympus HF Anschluss, Handgriff demontierbar, zusätzlicher Spülanschluss, autoklavierbar 134° C**

Polypectomy snare, half moon, multifilament wire, Olympus HF connection, removable handle, additional rinse connection, autoclavable 134° C

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Größe Ø snare	unsteril non-sterile	1 VE pack. unit
FE-2423DB20	2,4 x 2300 mm	2,8 mm	20 mm	o	2 Stk. + Griff / 2 pcs. + handle
FE-2423DB30	2,4 x 2300 mm	2,8 mm	30 mm	o	2 Stk. + Griff / 2 pcs. + handle

Hexagonale Schlingen auf Anfrage erhältlich! / Hexagonal snares available on request!