



**Steinfänger „Spirality“, 4-fach, Nitinol – nicht Lithotripsie geeignet, autoklavierbar 134° C**

Retrieval basket „Spirality Basket“, Nitinol – unsuitable for lithotripsy, autoclavable 134° C

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Korbgröße Ø basket	unsteril non-sterile	1 VE pack. unit
FE-2419GL220	2,4 x 1900 mm	2,8 mm	20 mm	o	1 Stk./pcs.
FE-2419GL230	2,4 x 1900 mm	2,8 mm	30 mm	o	1 Stk./pcs.
FE-2423GL230	2,4 x 2300 mm	2,8 mm	30 mm	o	1 Stk./pcs.



**Steinfänger „Ball Basket“, 4-fach, monofiler Draht – nicht Lithotripsie geeignet, autoklavierbar 134° C**

Retrieval basket „Ball Basket“, monofilament wire – unsuitable for lithotripsy, autoclavable 134° C

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Korbgröße Ø basket	unsteril non-sterile	1 VE pack. unit
FE-2419GL320	2,4 x 1900 mm	2,8 mm	20 mm	o	1 Stk./pcs.
FE-2423GL330	2,4 x 2300 mm	2,8 mm	30 mm	o	1 Stk./pcs.